


Anti-Phospho-JunD (S255) Rabbit Monoclonal Antibody For Wb And Elisa - P17535

Item Number: CM0001001



Introduction

Phospho-JunD (S255) Rabbit Monoclonal Antibody validated for WB and ELISA. Detects human and mouse phosphorylated JunD. Ideal for AP-1 signaling research. - P17535

[Learn More](#)

Parameter	Value
Product Name	Phospho-JunD-S255 Rabbit mAb
Remarks/Alias	AP-1; Phospho-JunD-S255
Species	Human
GeneID	3727
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S255 of human JunD (P17535).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P17535
Protein Weight	35kDa
Shipping	Ice bag