


Anti-Phospho-JunD-S255 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P17535

Item Number: CM0000464



Introduction

Phospho-JunD-S255 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Detects human and rat JunD phosphorylated at S255 (~35 kDa), a nuclear transcription factor of the AP-1 complex.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JunD-S255 Rabbit mAb
Remarks/Alias	AP-1; Phospho-JunD-S255
Species	Human
GeneID (Human)	3727
GeneID	3727
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S255 of human Phospho-JunD-S255 (NP_005345.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P17535
Protein Weight	35 kDa
Shipping	Ice bag