


Anti-Phospho-Kap1/Trim28-S824 Rabbit Monoclonal Antibody For Wb, Elisa - Q13263

Item Number: CM0001256



Introduction

Phospho-specific rabbit monoclonal antibody targeting human TRIM28 (KAP-1) phosphorylated at Ser-824 for WB and ELISA. Detects human and mouse TRIM28, a transcriptional repressor and E3 ligase involved in chromatin regulation and genomic stability. Ideal for signaling and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-KAP1/TRIM28-S824 Rabbit mAb
Remarks/Alias	KAP1; TF1B; RNF96; TIF1B; PPP1R157; TIF1beta; Phospho-KAP1/TRIM28-S824
Species	Human
Gene ID	10155
Immunogen	Synthetic phosphorylated peptide around S824 of human KAP1/TRIM28 (NP_005753.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13263
Protein Weight	79kDa/89kDa
Shipping	Ice bag