

Anti-Phospho-Kras-S89 Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - P01116

Item Number: CM0001852



Introduction

Phospho-KRAS-S89 Rabbit pAb detects KRAS phosphorylated at serine 89. Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for cancer signaling research and oncogene studies. One-stop sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Phospho-KRAS-S89 Rabbit pAb
Remarks/Alias	NS; NS3; OES; CFC2; RALD; K-Ras; KRAS1; KRAS2; RASK2; KI-RAS; C-K-RAS; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; K-Ras 2; 'C-K-RAS; c-Ki-ras; c-Ki-ras2; Phospho-KRAS-S89
Species	Human
Gene ID (Human)	3845
Gene ID	3845
Immunogen	A synthetic phosphorylated peptide around S89 of human KRAS (NP_001356715.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01116
Protein Weight	22kDa
Shipping	Ice bag