

Anti-Phospho-Limk2-T505 Rabbit Polyclonal Antibody For Wb And If/Icc - P53671

Item Number: CM0001824



Introduction

High specificity phospho-specific rabbit polyclonal antibody against LIMK2 at T505, validated for Western blot and immunofluorescence in human, mouse, and rat samples. Essential for actin dynamics and cell motility research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LIMK2-T505 Rabbit pAb
Remarks/Alias	LIMK2; Phospho-LIMK2-T505
Species	Human
GeneID (Human)	3985
GeneID	3985
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding T505 of human LIMK2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Mouse, Rat
SWISS	P53671
Protein Weight	72kDa
Shipping	Ice bag