


Anti-Phospho-Lrp6 (S1490) Rabbit Polyclonal Antibody For Wb, Elisa - O75581

Item Number: CM0001778



Introduction

Rabbit polyclonal antibody against phosphorylated LRP6 at Ser1490, for WB and ELISA applications. LRP6 is a Wnt co-receptor involved in beta-catenin signaling, crucial for development and disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LRP6-S1490 Rabbit pAb
Remarks/Alias	ADCAD2; STHAG7; Phospho-LRP6-S1490
Species	Human
Gene ID (Human)	4040
Gene ID	4040
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S1490 of human LRP6 (NP_002327.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75581
Protein Weight	180kDa
Shipping	Ice bag