

Anti-Phospho-Lrrk2-S935 Rabbit Monoclonal Antibody For Wb And Elisa - Q5S007

Item Number: CM0001970



Introduction

Rabbit monoclonal antibody detecting phosphorylated LRRK2 (Ser935) for WB and ELISA. Cross-reacts with mouse and rat. High specificity, validated with synthetic phosphopeptide. Ideal for Parkinson's disease and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-LRRK2-S935 Rabbit mAb
Remarks/Alias	PARK8; RPK7; ROCO2; AURA17; DARDARIN; Phospho-LRRK2-S935
Species	Human
Gene ID (Human)	120892
Immunogen	A synthetic phosphorylated peptide around S935 of human LRRK2 (Q5S007)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5S007
Protein Weight	286kDa
Shipping	Ice bag