

Anti-Phospho-Map2K3-S189/Map2K6-S207 Rabbit Polyclonal Antibody For Wb, Elisa - P46734 / P52564

Item Number: CM0001288

Introduction

Rabbit polyclonal antibody detecting phosphorylated MAP2K3 (S189) and MAP2K6 (S207) in human, mouse, and rat. Suitable for Western blot and ELISA. Ideal for studying stress-activated MAP kinase cascades and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K3-S189/MAP2K6-S207 Rabbit pAb
Remarks/Alias	MAP2K3/MAP2K6; Phospho-MAP2K3-S189/MAP2K6-S207
Species	Human
GeneID (human)	5606/5608
GeneID	5606/5608
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46734/P52564
Protein Weight	36kDa/39kDa/31kDa/37kDa
Shipping	Ice bag