


Anti-Phospho-Map2K4 (S257) Rabbit Monoclonal Antibody For Wb, Ihc-P - P45985

Item Number: CM0001357



Introduction

Phospho-MAP2K4 (S257) Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat phospho-MAP2K4, a key JNK/p38 MAPK kinase involved in stress signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K4-S257 Rabbit mAb
Remarks/Alias	JNKK; MEK4; MKK4; SEK1; SKK1; JNKK1; SERK1; MAPKK4; PRKMK4; SAPKK1; SAPKK-1
Species	Human
GeneID (Human)	6416
Immunogen	A synthetic phosphorylated peptide around S257 of human MAP2K4 (NP_003001.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45985
Protein Weight	44kDa
Shipping	Ice bag