

Anti-Phospho-Map3K5 (S1033) Rabbit Polyclonal Antibody For Wb And Elisa - Q99683

Item Number: CM0000958



Introduction

Rabbit polyclonal antibody specific to phosphorylated MAP3K5 (ASK1) at Ser1033. Validated for WB and ELISA in human samples. Detects stress-related kinase involved in apoptosis and MAP kinase signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP3K5-S1033 Rabbit pAb
Remarks/Alias	ASK1; MEKK5; MAPKKK5; Phospho-MAP3K5-S1033
Species	Human
GeneID (Human)	4217
GeneID	4217
Immunogen	A synthetic phosphorylated peptide around S1033 of human MAP3K5 (NP_005914.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99683
Protein Weight	155kDa
Shipping	Ice bag