

Anti-Phospho-Map3K5-T918 Rabbit Polyclonal Antibody For Wb, Elisa - Q99683

Item Number: CM0001401



Introduction

Rabbit polyclonal antibody targeting phosphorylated Thr918 of human MAP3K5 (ASK1, MEKK5). Validated for Western Blot and ELISA. Ideal for studying apoptosis, oxidative stress response, and MAP kinase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP3K5-T918 Rabbit pAb
Remarks/Alias	ASK1; MEKK5; MAPKKK5; Phospho-MAP3K5-T918
Species	Human
GeneID	4217
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99683
Protein Weight	155kDa
Shipping	Ice bag