


Anti-Phospho-Mapk14 (T180/Y182) / Mapk11 (T180/Y182) / Mapk13 (T180/Y182) / Mapk12 (T183/Y1845) Rabbit Polyclonal Antibody For Wb, Elisa - Q15759/P53778/O15264/Q16539

Item Number: CM0001706



Introduction

This rabbit polyclonal antibody detects phosphorylated p38 MAPKs (MAPK11, MAPK12, MAPK13, MAPK14) at key activation sites, with WB and ELISA applications and human, rat cross-reactivity, ideal for studying stress signaling and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAPK14-T180/Y182+MAPK11-T180/Y182+MAPK13-T180/Y182+MAPK12-T183/Y1845 Rabbit pAb
Species	Human
GenelD (Human)	5600/6300/5603/1432
GenelD	5600/6300/5603/1432
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T180 & Y182 of human MAPK14MAPK11 (NP_002742.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15759/P53778/O15264/Q16539
Shipping	Ice bag