

Anti-Phospho-Mapkapk-2/Mk2-T334 Rabbit Polyclonal Antibody For Ihc-P, Elisa - P49137

Item Number: CM0000836



Introduction

Phospho-MAPKAPK-2/MK2-T334 rabbit polyclonal antibody validated for IHC-P and ELISA, targeting human, mouse, and rat MAPKAPK2 phosphorylated at threonine 334. Ideal for stress signaling and cytokine regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAPKAPK-2/MK2-T334 Rabbit pAb
Remarks/Alias	MK2; MK-2; MAPKAP-K2; Phospho-MAPKAPK-2/MK2-T334
Species	Human
GeneID (Human)	9261
GeneID	9261
Immunogen	A synthetic phosphorylated peptide around T334 of human MAPKAPK2 (NP_116584.2?).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49137
Protein Weight	46kDa
Shipping	Ice bag