

Anti-Phospho-Mapkapk2 (T334) Rabbit Polyclonal Antibody For Wb, Elisa - P49137

Item Number: CM0001195



Introduction

Phospho-specific rabbit polyclonal antibody against MAPKAPK2 at threonine 334. Suitable for Western blot and ELISA; reacts with human, mouse, and rat MAPKAPK2. Ideal for stress signaling, inflammation, and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAPKAPK-2/MK2-T334 Rabbit pAb
Remarks/Alias	MK2; MK-2; MAPKAP-K2; Phospho-MAPKAPK-2/MK2-T334
Species	Human
Gene ID (Human)	9261
Gene ID	9261
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T334 of human MAPKAPK2 (NP_116584.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49137
Protein Weight	46kDa
Shipping	Ice bag