

Anti-Phospho-Marcks-S158 Rabbit Polyclonal Antibody For Wb - P29966

Item Number: CM0001233

Introduction

Rabbit polyclonal antibody specific for MARCKS phosphorylated at Ser158. Validated for Western blot detection in human, mouse, and rat. Ideal for studying PKC signaling, actin cytoskeleton dynamics, and inflammation.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-MARCKS-S158 Rabbit pAb
Remarks/Alias	MACS; 80K-L; PKCSL; PRKCSL; Phospho-MARCKS-S158
Species	Human
GeneID (Human)	4082
GeneID	4082
Immunogen	Synthetic peptide corresponding to residues surrounding S158 of human MARCKS (phospho-specific)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P29966
Protein Weight	32 kDa
Shipping	Ice bag