

Anti-Phospho-Mek1-S298 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q02750

Item Number: CM0000967



Introduction

Rabbit polyclonal antibody against phosphorylated MEK1 (S298), validated for WB, IP, ELISA in human, mouse, rat. Ideal for MAPK/ERK pathway studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK1-S298 Rabbit pAb
Remarks/Alias	MEL; CFC3; MEK1; MKK1; MAPKK1; PRKMK1; Phospho-MEK1-S298
Species	Human
Gene ID	5604
Immunogen	A synthetic phosphorylated peptide around S298 of human MEK1 (NP_002746.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02750
Protein Weight	43kDa
Shipping	Ice bag