

Anti-Phospho-Mek2 (T394) Rabbit Polyclonal Antibody For Wb And If/Icc - P36507

Item Number: CM0001573



Introduction

This rabbit polyclonal antibody detects endogenous levels of MEK2 only when phosphorylated at Threonine 394, a critical kinase in the MAPK/ERK signaling cascade. Suitable for Western blot and immunofluorescence applications in human cells.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MEK2-T394 Rabbit pAb
Remarks/Alias	CFC4; MEK2; MKK2; MAPKK2; PRKMK2; Phospho-MEK2-T394
Species	Human
Gene ID (Human)	5605
Gene ID	5605
Immunogen	A phospho specific peptide corresponding to residues surrounding T394 of human MEK2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human
SWISS	P36507
Protein Weight	44kDa
Shipping	Ice bag