

Anti-Phospho-Met-Y1234/1235 Rabbit Pab For Wb, Elisa - P08581

Item Number: CM0001305



Introduction

Anti-Phospho-MET-Y1234/1235 Rabbit pAb from CamelBio detects phosphorylated MET (Y1234/1235) in human samples via WB and ELISA. Ideal for studying HGF/MET signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MET-Y1234/1235 Rabbit pAb
Remarks/Alias	DA11; HGFR; AUTS9; RCCP2; c-Met; DFNB97; Phospho-MET-Y1234/1235
Species	Human
GeneID	4233 (human)
Immunogen	A synthetic phosphorylated peptide around Y1234 & Y1235 of human MET (NP_000236.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08581
Protein Weight	156kDa
Shipping	Ice bag