

Anti-Phospho-Mkk3 (S189) + Mkk6 (S207) Rabbit Monoclonal Antibody For Wb/If/Elisa - P46734 / P52564

Item Number: CM0000820

Introduction

Phospho-specific rabbit monoclonal antibody targeting phosphorylation at S189 of MKK3 and S207 of MKK6. Validated for Western blot, IF/ICC, and ELISA in human and mouse samples to study stress-activated MAPK pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MKK3-S189+MKK6-S207 Rabbit mAb
Remarks/Alias	MAP2K3/MAP2K6; Phospho-MKK3-S189+MKK6-S207
Species	Human
Gene ID	5606/5608
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S189/S207 of human MKK3/MKK6.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46734/P52564
Protein Weight	36kDa/39kDa/31kDa/37kDa
Shipping	Ice bag