

Anti-Phospho-Mst1-T183+Mst2-T180 Rabbit Monoclonal Antibody For Wb And Elisa - Q13188 / Q13043

Item Number: CM0001324



Introduction

Phospho-specific rabbit monoclonal antibody recognizing MST1 (T183) and MST2 (T180) in human, mouse, rat. Ideal for WB and ELISA in Hippo pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MST1-T183+MST2-T180 Rabbit mAb
Remarks/Alias	KRS1; MST2; KRS2; MST1
Species	Human
GeneID (human)	6788/6789
GeneID	6788/6789
Immunogen	A synthetic phosphorylated peptide around T180/T183 of human MST2 (NP_006272.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13188/Q13043
Protein Weight	56kDa
Shipping	Ice bag