


Anti-Phospho-Mst1/Stk4-T387 Rabbit Polyclonal Antibody For Western Blot And Elisa - Q13043

Item Number: CM0001244



Introduction

Rabbit polyclonal antibody against phosphorylated MST1/STK4 at T387, validated for Western blot and ELISA with cross-reactivity to human, mouse, and rat samples. Suitable for studying apoptosis and the Hippo signaling pathway.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MST1/STK4-T387 Rabbit pAb
Remarks/Alias	KRS2; MST1; YSK3; Phospho-MST1/STK4-T387
Species	Human
Gene ID (Human)	6789
Immunogen	A synthetic phosphorylated peptide around T387 of human MST1/STK4 (NP_006273.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13043
Protein Weight	56kDa
Shipping	Ice bag