

Anti-Phospho-Mtor-S2448 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P42345

Item Number: CM0001894



Introduction

Phospho-mTOR (Ser2448) rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat mTOR, targeting a key phosphorylation site in mTOR signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-mTOR-S2448 Rabbit mAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; Phospho-mTOR-S2448
Species	Human
GeneID (Human)	2475
GeneID	2475
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag