

Anti-Phospho-Mtor-S2448 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42345

Item Number: CM0000987



Introduction

Phospho-mTOR-S2448 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat mTOR when phosphorylated at serine 2448. Ideal for mTOR signaling pathway analysis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-mTOR-S2448 Rabbit pAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; Phospho-mTOR-S2448
Species	Human
GeneID	2475
Immunogen	A synthetic phosphorylated peptide around S2448 of human mTOR (NP_004949.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag