

Anti-Phospho-Mtor (S2481) Rabbit Polyclonal Antibody For Wb, Elisa - P42345

Item Number: CM0000582



Introduction

Validate mTOR signaling with this rabbit polyclonal antibody specific to phosphorylated Ser2481. Suitable for WB and ELISA in human, mouse, and rat samples. Ideal for studying metabolic and growth pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-mTOR-S2481 Rabbit pAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; Phospho-mTOR-S2481
Species	Human
Gene ID (Human)	2475
Gene ID	2475
Immunogen	A synthetic phosphorylated peptide around S2481 of human mTOR (NP_004949.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag