

Anti-Phospho-Ncf4 (P40-Phox) T154 Rabbit Polyclonal Antibody For Wb And Elisa - Q15080

Item Number: CM0001881



Introduction

Phospho-specific rabbit polyclonal antibody against human NCF4/p40-phox (T154 phosphorylation site). Suitable for WB and ELISA, it aids in investigating NADPH oxidase complex activation and phagocyte biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NCF4/p40-phox-T154 Rabbit pAb
Remarks/Alias	NCF; CGD3; P40PHOX; SH3PXD4; Phospho-NCF4/p40-phox-T154
Species	Human
GeneID (Human)	4689
GeneID	4689
Immunogen	A synthetic phosphorylated peptide around T154 of human NCF4 (NP_000622.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15080
Protein Weight	39kDa
Shipping	Ice bag