


Anti-Phospho-NdrG1-S330 Rabbit Polyclonal Antibody For Wb, Elisa - Q92597

Item Number: CM0001476



Introduction

Phospho-NDRG1 (S330) rabbit polyclonal antibody validated for Western blot (WB) and ELISA, reactive with human and mouse. Useful for stress response, tumor suppression, and myelination research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NDRG1-S330 Rabbit pAb
Remarks/Alias	GC4; RTP; DRG1; NDR1; NMSL; TDD5; CAP43; CMT4D; DRG-1; HMSNL; RIT42; TARG1; PROXY1; Phospho-NDRG1-S330
Species	Human
GeneID (Human)	10397
GeneID	10397
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S330 of human NDRG1 (NP_001128714.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92597
Protein Weight	43kDa
Shipping	Ice bag