


Anti-Phospho-NdrG1 (T346) Rabbit Polyclonal Antibody For Wb And Elisa - Q92597

Item Number: CM0001878



Introduction

Rabbit polyclonal antibody against Phospho-NDRG1 (T346) for WB and ELISA applications. Detects the stress-responsive tumor suppressor NDRG1 phosphorylated at threonine 346. Validated in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NDRG1-T346 Rabbit pAb
Remarks/Alias	GC4, RTP, DRG1, NDR1, NMSL, TDD5, CAP43, CMT4D, DRG-1, HMSNL, RIT42, TARG1, PROXY1, Phospho-NDRG1-T346
Species	Human
Gene ID	10397
Immunogen	Synthetic phosphorylated peptide around T346 of human NDRG1 (NP_001128714.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92597
Protein Weight	43kDa
Shipping	Ice bag