


Anti-Phospho-Nf-Kb P65/Rela (S536) Monoclonal Antibody For Wb, If/Icc, Elisa - Q04206

Item Number: CM0001427



Introduction

Phospho-specific rabbit monoclonal antibody against human NF-κB p65 phosphorylated at Ser536. Cross-reacts with mouse and rat. For WB, IF/ICC, ELISA. Ideal for studying NF-κB activation in inflammation and immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB p65/RelA-S536 Rabbit mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-κB p65/RelA-S536
Species	Human
Gene ID (Human)	5970
Gene ID	5970
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S536 of human NF-κB (NP_068810.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	65kDa
Shipping	Ice bag