


Anti-Phospho-Nf-Kb P65/Rela (Ser529) Rabbit Mab For Wb, Elisa - Q04206

Item Number: CM0001248



Introduction

Phospho-specific rabbit monoclonal antibody targeting the S529 phosphorylation site of human NF- κ B p65 (RelA), validated for Western blot (WB) and ELISA. Ideal for studying NF- κ B signaling in inflammation, immunity, and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF- κ B p65/RelA-S529 Rabbit mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF- κ B p65/RelA-S529
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S529 of human NF- κ B p65/RelA (NP_068810.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag