


Anti-Phospho-Nf-Kb P65/Rela-T435 Rabbit Polyclonal Antibody For Wb And Elisa - Q04206

Item Number: CM0001187



Introduction

Rabbit polyclonal antibody against phosphorylated NF-kB p65 (RelA) at Thr435. Validated for WB and ELISA, reacts with human and mouse samples. Ideal for NF-kB signaling and inflammation research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-NF-kB p65/RelA-T435 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-T435
Species	Human
Gene ID (Human)	5970
Gene ID	5970
Immunogen	Synthetic phosphorylated peptide around T435 of human RelA (NP_068810.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag