


Anti-Phospho-Nf-Kb2-S866 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q00653

Item Number: CM0001456



Introduction

Phospho-NF-κB2-S866 Rabbit Polyclonal Antibody specifically detects endogenous phosphorylated NF-κB2 p100/p52 at Ser866. Ideal for WB, IHC-P, IF/ICC, and ELISA applications in human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB2-S866 Rabbit pAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; Phospho-NF-κB2-S866
Species	Human
Gene ID (Human)	4791
Immunogen	A synthetic phosphorylated peptide around S866 of human NF-κB2 (NP_001070962.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag