


Anti-Phospho-Nfkb1-S337 Rabbit Mab For Wb, Elisa - P19838

Item Number: CM0000738



Introduction

Phospho-NFKB1-S337 Rabbit Monoclonal Antibody detects endogenous levels of NFKB1 only when phosphorylated at serine 337. Validated for Western blot and ELISA in human samples—ideal for inflammation and immunity research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-NFKB1-S337 Rabbit mAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta
Species	Human
Gene ID (Human)	4790
Gene ID	4790
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S337 of human NFKB1 (NP_001158884.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19838
Protein Weight	85kDa/105kDa
Shipping	Ice bag