

Anti-Phospho-Nfkb1-S893 Rabbit Polyclonal Antibody For If/Icc - P19838

Item Number: CM0000944



Introduction

Phospho-specific rabbit polyclonal antibody targeting NFKB1 p105 phosphorylated at S893. Validated for immunofluorescence/immunocytochemistry in human samples. Supports research on NF-κB signaling, inflammation, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NFKB1-S893 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-κB; CVID12; NF-κB1; NFKB-p50; NFKappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NFKB1-S893
Species	Human
GeneID (Human)	4790
GeneID	4790
Immunogen	Synthetic peptide corresponding to residues surrounding S893 of human NFKB1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IF/ICC
Cross Reactivity	Human
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag