


Anti-Phospho-Nfkb1 (S932) Rabbit Polyclonal Antibody For Wb - P19838

Item Number: CM0000536



Introduction

Rabbit polyclonal antibody specifically recognizing NFKB1 phosphorylated at serine 932. Validated for Western blot in human, mouse, and rat samples. Ideal for NF-kB signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NFKB1-S932 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFKappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NFKB1-S932
Species	Human
GeneID (Human)	4790
GeneID	4790
Immunogen	Synthetic peptide; A phospho specific peptide corresponding to residues surrounding S932 of human NFKB1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag