

Anti-Phospho-Nmdar1-S890 Rabbit Polyclonal Antibody For Wb, Elisa - Q05586

Item Number: CM0000533



Introduction

Phospho-specific rabbit polyclonal antibody detecting NMDAR1 phosphorylated at Ser890. Validated in WB and ELISA, with cross-reactivity to human and rat. Ideal for studying NMDA receptor signaling in synaptic plasticity and neuropathology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NMDAR1-S890 Rabbit pAb
Remarks/Alias	NR1; MRD8; GluN1; NMDA1; DEE101; NDHMSD; NDHMSR; NMD-R1; NMDAR1; Phospho-NMDAR1-S890
Species	Human
Gene ID (Human)	2902
Gene ID	2902
Immunogen	A synthetic phosphorylated peptide around S890 of human NMDAR1 (NP_015566.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q05586
Protein Weight	105kDa
Shipping	Ice bag