

Anti-Phospho-Nmdar1 (S896/S897) Polyclonal Antibody For Wb, Elisa - Q05586

Item Number: CM0001674



Introduction

Rabbit polyclonal antibody against phosphorylated Ser896/Ser897 of human NMDAR1 (GluN1). Suitable for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for synaptic plasticity and NMDA receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NMDAR1-S896/S897 Rabbit pAb
Remarks/Alias	NR1; MRD8; GluN1; NMDA1; DEE101; NDHMSD; NDHMSR; NMD-R1; NMDAR1; Phospho-NMDAR1-S896/S897
Species	Human
GeneID (Human)	2902
GeneID	2902
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S896 & S897 of human NMDAR1 (NP_015566.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05586
Protein Weight	105kDa
Shipping	Ice bag