


Anti-Phospho-Npm1-T199 Rabbit Polyclonal Antibody For Wb, Elisa - P06748

Item Number: CM0001752



Introduction

Rabbit polyclonal antibody detecting phosphorylated NPM1 at Thr199 in human samples. Suitable for Western blot and ELISA. Ideal for studying nucleophosmin phosphorylation in ribosome biogenesis, centrosome duplication, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NPM1-T199 Rabbit pAb
Remarks/Alias	B23; NPM; Phospho-NPM1-T199
Species	Human
GeneID (Human)	4869
Immunogen	Synthetic phosphorylated peptide around T199 of human NPM1 (NP_001032827.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06748
Protein Weight	33kDa
Shipping	Ice bag