


Anti-Phospho-Nrf2-S40 Rabbit Monoclonal Antibody For Wb - Q16236

Item Number: CM0001974



Introduction

Phospho-NRF2-S40 Rabbit mAb detects phosphorylated NRF2 in human, mouse, and rat. Suitable for WB, IF, IHC-P, and ELISA. Ideal for oxidative stress and cytoprotective gene research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NRF2-S40 Rabbit mAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH
Species	Human
GeneID	4780
Immunogen	A synthetic phosphorylated peptide around S40 of human NRF2 (Q16236)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag