


Anti-Phospho-Nrf2-S40 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16236

Item Number: CM0000604



Introduction

Validated for WB, IF/ICC, ELISA, this rabbit polyclonal antibody detects human NRF2 phosphorylated at serine 40. An essential tool for oxidative stress, redox signaling, and cytoprotection research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NRF2-S40 Rabbit pAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH; Phospho-NRF2-S40
Species	Human
Gene ID (Human)	4780
Gene ID	4780
Immunogen	A synthetic phosphorylated peptide around S40 of human NRF2 (NP_006155.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q16236
Protein Weight	68kDa
Shipping	Ice bag