

Anti-Phospho-Nucleolin-T76 Rabbit Monoclonal Antibody For IHC-P, Elisa - P19338

Item Number: CM0001614



Introduction

Phospho-Nucleolin-T76 Rabbit Monoclonal Antibody for detecting human Nucleolin phosphorylated at T76. Validated for IHC-P and ELISA. Supports research on ribosome biogenesis and chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Nucleolin-T76 Rabbit mAb
Remarks/Alias	C23; Nsr1; Phospho-Nucleolin-T76
Species	Human
GeneID (Human)	4691
GeneID	4691
Immunogen	Synthetic peptide corresponding to residues surrounding T76 of Human Nucleolin.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P19338
Protein Weight	77kDa
Shipping	Ice bag