


Anti-Phospho-P38 Mapk (T180) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16539 / P53778 / Q15759

Item Number: CM0001036



Introduction

Rabbit monoclonal antibody targeting phosphorylated p38 MAPK (Thr-180), validated for WB, IHC-P, ELISA with reactivity across human, mouse, and rat. Essential for studying stress-activated MAP kinase signaling, inflammation, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-T180 Rabbit mAb
Remarks/Alias	MAPK14; CSBP; CSBP1; CSBP2; CSPB1; EXIP; Mxi2; PRKM14; PRKM15; RK; SAPK2A; p38; p38ALPHA; p38 MAPK; Phospho-p38 MAPK-T180
Species	Human
GeneID (human)	1432, 6300, 5600
GeneID	1432
Immunogen	A synthetic phosphorylated peptide around T180 of human p38 MAPK (NP_620581.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539, P53778, Q15759
Protein Weight	41kDa
Shipping	Ice bag