


Anti-Phospho-P38 Mapk (T180/Y182) Rabbit Monoclonal Antibody For Wb And Elisa - Q15759 / P53778 / O15264 / Q16539

Item Number: CM0000843



Introduction

Phospho-specific rabbit monoclonal antibody against phosphorylated p38 MAPK (T180/Y182). Validated for WB and ELISA, cross-reactive with human, mouse, and rat. Ideal for stress and inflammatory signaling pathway research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-p38 MAPK-T180/Y182 Rabbit mAb
Remarks/Alias	MAPK14; CSBP; CSBP1; CSBP2; CSPB1; EXIP; Mxi2; PRKM14; PRKM15; RK; SAPK2A; p38; p38ALPHA; p38 MAPK; Phospho-p38 MAPK-T180/Y182
Species	Human
GenelD (Human)	5600/6300/5603/1432
GenelD	5600/6300/5603/1432
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T180 & Y182 of human MAPK (NP_002742.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15759/P53778/O15264/Q16539
Protein Weight	29kDa/34kDa/35kDa/41kDa
Shipping	Ice bag