


Anti-Phospho-P38 Mapk (Y182) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q16539

Item Number: CM0001072



Introduction

Phospho-specific rabbit monoclonal antibody against p38 MAPK (MAPK14) phosphorylated at tyrosine 182 (Y182). Validated for WB, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for studying stress and inflammatory signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-Y182 Rabbit mAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y182 of human p38 MAPK (NP_620581.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag