

Anti-Phospho-P38 Mapk (Y182) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16539

Item Number: CM0002001



Introduction

Phospho-specific rabbit monoclonal antibody targeting p38 MAPK phosphorylated at Y182; validated for WB, IHC-P, IF/ICC, ELISA; detects human, mouse, rat MAPK14 (Q16539), key in stress, inflammation, and apoptosis signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-Y182 Rabbit mAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; Phospho-p38 MAPK-Y182
Species	Human
Gene ID (Human)	1432
Immunogen	A synthetic phosphorylated peptide around Y182 of human MAPK
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag