

Anti-Phospho-P38 Mapk (Y182) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16539

Item Number: CM0000713



Introduction

High-specificity rabbit polyclonal antibody targeting phosphorylated p38 MAPK (Tyr182). Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for investigating stress and inflammatory signal transduction pathways.

[Learn More](#)

Parameter	Description
Product Name	Phospho-p38 MAPK-Y182 Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; Phospho-p38 MAPK-Y182
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y182 of human MAPK (NP_620581.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag