

## Anti-Phospho-P53-S15 Polyclonal Antibody For Wb, If/Icc, Elisa - P04637

Item Number: CM0000917



### Introduction

Rabbit polyclonal antibody specific for phosphorylated p53 at Ser15. Validated for Western blot, IF/ICC, and ELISA. For studying DNA damage response, apoptosis, and cell cycle regulation in human, mouse, rat samples.

[Learn More](#)

| Parameter        | Value                                                                                       |
|------------------|---------------------------------------------------------------------------------------------|
| Product Name     | Phospho-p53-S15 Rabbit pAb                                                                  |
| Remarks/Alias    | P53; BCC7; LFS1; BMF55; TRP53; Phospho-p53-S15                                              |
| Species          | Human                                                                                       |
| GeneID (human)   | 7157                                                                                        |
| GeneID           | 7157                                                                                        |
| Immunogen        | Synthetic peptide: A synthetic phosphorylated peptide around S15 of human p53 (NP_000537.3) |
| Source           | Rabbit                                                                                      |
| Category         | Phosphorylated Antibodies                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                           |
| SWISS            | P04637                                                                                      |
| Protein Weight   | 44kDa                                                                                       |
| Shipping         | Ice bag                                                                                     |