

Anti-Phospho-P53-S315 Polyclonal Antibody For Wb, Elisa - P04637

Item Number: CM0001375



Introduction

Rabbit polyclonal antibody specifically recognizing phosphorylated p53 at Ser315. Validated for WB and ELISA in human and rat samples. Suitable for apoptosis, cell cycle, and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S315 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMF55; TRP53; Phospho-p53-S315
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S315 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag