


Anti-Phospho-P53 (S33) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P04637

Item Number: CM0001860



Introduction

Phospho-p53-S33 Rabbit pAb detects endogenous levels of p53 phosphorylated at Ser33 in human, mouse, and rat. Suitable for WB, IHC-P, IP, and ELISA, it is a valuable tool for studying p53-mediated cell cycle arrest, apoptosis, and tumor suppression signaling.

[Learn More](#)

Parameter	Details
Product Name	Phospho-p53-S33 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S33
Species	Human
Gene ID (Human)	7157
Gene ID	7157
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around S33 of human p53 (NP_000537.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag