


Anti-Phospho-P53-S392 Rabbit Polyclonal Antibody For Wb, Elisa - P04637

Item Number: CM0000927



Introduction

Phospho-p53-S392 rabbit polyclonal antibody (P04637) validated for WB and ELISA applications. Suitable for human, mouse, rat samples. Ideal for studying p53 signaling, apoptosis, and tumor suppression in cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S392 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S392
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	Synthetic Phosphopeptide: A synthetic phosphorylated peptide around S392 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag