


Anti-Phospho-P53 (S9) Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P04637

Item Number: CM0001773



Introduction

Phospho-p53 (S9) Rabbit Polyclonal Antibody (pAb) for Western blotting (WB), immunoprecipitation (IP), and ELISA. Recognizes human p53 phosphorylated at serine 9. Suitable for research on tumor suppression, apoptosis, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S9 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S9
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	A synthetic phosphorylated peptide around S9 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag